



HBR



SMS



BA9

BIOFREEDOMTM
POLYMER- & CARRIER-FREE DRUG-COATED CORONARY STENT SYSTEM

Powered by



Freedom to treat

LEADERS **FREE**

**The Landmark trial evaluating
BioFreedom in High Bleeding
Risk (HBR) patients with
1 month DAPT**

Polymer-free Drug-Coated Coronary Stents
in Patients at High Bleeding Risk



BIOSENSORS
INTERNATIONALTM

The Landmark trial evaluating BioFreedom in High Bleeding Risk (HBR) patients with 1 month DAPT

LEADERS FREE Trial Design

Prospective, double-blind randomized (1:1) trial
2466 High Bleeding Risk (HBR) PCI patients

BioFreedomTM
DCS

vs.

GazelleTM
BMS

DAPT mandated for 1 month only, followed by long-term SAPT

Primary safety endpoint

Composite of cardiac death, MI, definite / probable stent thrombosis at 1 year (non-inferiority then superiority)

Primary efficacy endpoint

Clinically-driven TLR at 1 year (superiority)

A Significantly Sicker Population

| Baseline Characteristics | DCS (%) | BMS (%) |
|---------------------------------------|------------|------------|
| Mean age | 75.7 + 9.4 | 75.7 + 9.3 |
| Female gender | 29.8 | 30.9 |
| BMI | 27.5 ± 4.8 | 27.2 ± 4.6 |
| Diabetes | 34.0 | 32.3 |
| NSTEMI presentation | 22.4 | 23.2 |
| STEMI presentation | 4.7 | 4.0 |
| Prior MI | 19.6 | 21.4 |
| Prior PCI | 22.2 | 21.9 |
| Prior CABG | 9.4 | 10.1 |
| Multivessel CAD | 62.9 | 61.6 |
| Congestive heart failure | 14.4 | 12.4 |
| Atrial fibrillation | 34.9 | 34.6 |
| Peripheral vascular disease | 15.7 | 15.8 |
| Chronic obstructive pulmonary disease | 10.9 | 11.7 |

No significant differences in baseline characteristics between the groups


- Average age is 76 years (10 years older than typical all-comers trial population)
- 62% multivessel disease
- 35% atrial fibrillation
- 13% congestive heart failure
- 33% diabetes mellitus

Values reflect the mean of the aggregate population

Defining the HBR population

Neglected and Systematically Excluded from Stent Trials


On average, patients in LEADERS FREE had 1.7 inclusion criteria



Significantly Safer than BMS

29% Reduction in the Rate of the Composite of Cardiac Death, MI, ST

Primary Safety Endpoint (Composite of Cardiac Death, MI, ST)



Number at Risk

| | 0 | 90 | 180 | 270 | 390 |
|------------|------|------|------|------|------|
| DCS | 1221 | 1146 | 1105 | 1081 | 1045 |
| BMS | 1211 | 1115 | 1066 | 1037 | 1000 |


390 days chosen for assessing primary endpoint to capture potential events driven by the 360 day follow up contact
Hazard Ratio (HR) 0.71; 95% CI : 0.56 to 0.91; P = 0.005 for superiority

* Relative Risk Reduction

Significantly more Effective than BMS

50% Reduction in the Rate of Restenosis.

Primary Efficacy Endpoint (Clinically-Driven TLR)



Number at Risk

| | 0 | 90 | 180 | 270 | 390 |
|------------|------|------|------|------|------|
| DCS | 1221 | 1167 | 1130 | 1098 | 1053 |
| BMS | 1211 | 1131 | 1072 | 1034 | 984 |

390 days chosen for assessing primary endpoint to capture potential events driven by the 360 day follow up contact
Hazard Ratio (HR) 0.50; 95% CI : 0.37 to 0.69; P < 0.001 for superiority

* Relative Risk Reduction



HBR



SMS



BA9

BIOFREEDOMTM
POLYMER- & CARRIER-FREE DRUG-COATED CORONARY STENT SYSTEM

Powered by



Freedom to treat

BioFreedom is the only active stent with 1 month DAPT that has demonstrated superior outcomes to BMS



Significantly Safer than BMS



Significantly more Effective than BMS

With LEADERS FREE, BioFreedom becomes the **standard of care for High Bleeding Risk (HBR) patients**

Reference:

Urban et al. New England Journal of Medicine 2015; published ahead of print October 14. DOI: 10.1056/NEJMoa1503943

BioFreedomTM drug-coated stent and GazelleTM coronary stent are CE Mark approved. Data on file at Biosensors International for any sustained claims in this brochure.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

BioFreedom, Gazelle, Biolimus A9 and BA9 are trademarks or registered trademarks of Biosensors International Group, Ltd. All cited trademarks are the property of their respective owners

Not available for sale in the United States and certain other countries.

© 2015. Biosensors International Group, Ltd. All rights reserved.

www.biosensors.com



BIOSENSORS EUROPE SA

Rue de Lausanne 29
1110 Morges
Switzerland
Tel: +41 (0)21 804 80 00
Fax: +41 (0)21 804 80 01

BIOSENSORS INTERVENTIONAL TECHNOLOGIES PTE LTD

36 Jalan Tukang
Singapore 619 266
Tel: +65 6213 5777
Fax: +65 6213 5737