EN ENGLISH – INSTRUCTIONS FOR USE	7
FR FRANÇAIS – INSTRUCTION D'UTILISATION	9
DE DEUTSCH – GEBRAUCHSANWEISUNG	11
IT ITALIANO – ISTRUZIONI PER L'USO	14
ES ESPAÑOL – INSTRUCCIONES DE USO	16
PT PORTUGUÊS – INSTRUÇÕES DE UTILIZAÇÃO	18
NL NEDERLANDS – GEBRUIKSAANWIJZING	20
DA DANSK – BRUGSANVISNING	23
FI SUOMI – KÄYTTÖOHJEET	25
SV SVENSKA – BRUKSANVISNING	27
NO NORSK – BRUKSANVISNING	29
EL ΕΛΛΗΝΙΚΑ – ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	31
BG БЪЛГАРСКИ – ИНСТРУКЦИИ ЗА УПОТРЕБА	34
ET EESTI KEEL – KASUTUSJUHEND	37
LT LIETUVIŲ KALBA – NAUDOJIMO INSTRUKCIJA	39
LV LATVIEŠU VALODA – LIETOŠANAS INSTRUKCIJA	41
PL POLSKI – INSTRUKCJA UŻYCIA	43
CS ČEŠTINA – NÁVOD K POUŽITÍ	46
SK SLOVENČINA – NÁVOD NA POUŽITIE	48
HU MAGYAR – HASZNÁLATI UTASÍTÁS	50
RO ROMÂNĂ – INSTRUCȚIUNI DE UTILIZARE	52
HR HRVATSKI – UPUTE ZA UPOTREBU	54
SL SLOVENŠČINA – NAVODILA ZA UPORABO	56
TR TÜRKÇE – KULLANMA TALİMATI	58
MK МАКЕДОНСКИ – УПАТСТВА ЗА УПОТРЕБА	60
SR SRPSKI – UPUTSTVO ZA UPOTREBU	63
RU РУССКИЙ – ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ	65
ID BAHASA INDONESIA – PETUNJUK PENGGUNAAN	68
VI TIẾNG VIỆT – HƯỚNG DẪN SỬ DỤNG	70
ZH 简体中文 – 使用说明书	72
ZH 繁體中文 – 使用說明	74

EN SYMBOLS USED IN LABELING / **FR** SYMBOLES UTILISÉS POUR L'ÉTIQUETAGE / **DE** ZUR KENNZEICHNUNG VERWENDETE SYMBOLE / **IT** SIMBOLI UTILIZZATI PER L'ETICHETTATURA / **ES** SÍMBOLOS UTILIZADOS EN LAS ETIQUETAS / **PT** SÍMBOLOS UTILIZADOS NO RÓTULO / **NL** OP HET ETIKET GEBRUIKTE SYMBOLEN / **DA** ANVENDTE SYMBOLER PÅ ETIKETTER / **FI** MERKINNÖISSÄ KÄYTETYT SYMBOLIT / **SV** SYMBOLERI MÄRKNINGEN / **NO** SYMBOLER SOM BRUKES VED MERKING / **EL** ΣΥΜΒΟΛΑ ΠΟΥ ΧΡΗΣΙΜΟΠΟΙΟΥΝΤΑΙ ΣΤΙΣ ΕΤΙΚΕΤΕΣ / **BG** СИМВОЛИ ВЪРХУ ЕТИКЕТИТЕ / **ET** SILDIDEL KASUTATUD SÜMBOLID / **LT** ETIKETĖJE NAUDOJAMI SIMBOLIAI / **LV** MARKĒJUMĀ LIETOTIE SIMBOLI / **PL** SYMBOLE UŻYTE NA ETYKIECIE / **CS** SYMBOLY POUŽITÉ NA ŠTÍTKU / **SK** SYMBOLY POUŽITÉ NA ŠTÍTKU / **HU** A CÍMKÉN SZEREPLŐ SZIMBÓLUMOK / **RO** SIMBOLURI UTILIZATE LA ETICHETARE / **HR** SIMBOLIKOJISE UPOTREBLJAVAJUZA OZNAČIVANJE / **SL** SIMBOLI, UPORABLJENINA OZNAKAH / **TR** ETİKETLEMEDE KULLANILAN SEMBOLLER / **MK** СИМБОЛИ КОИ СЕ УПОТРЕБУВААТ ПРИ ЕТИКЕТИРАЊЕТО / **SR** SIMBOLI KOJI SE KORISTE NA NALEPNICAMA / **RU** УСЛОВНЫЕ ОБОЗНАЧЕНИЯ НА ЭТИКЕТКАХ / **ID** SIMBOL YANG DIGUNAKAN DALAM PELABELAN / **VI** KÝ HIỆU DÙNG GHI NHÃN / **ZH-CN** 标签所用符号 / **ZH-TW** 標籤所用的符號

	EN Manufacturer / FR Fabricant / DE Hersteller / IT Produttore / ES Fabricante / PT Fabricante / NL Fabrikant / DA Producent / FI Valmistaja / SV Tillverkare / NO Produsent / EL Κατασκευαστής / BG Производител / ET Tootja / LT Gamintojas / LV Ražotājs / PL Producent / CS Výrobce / SK Výrobca / HU Gyártó / RO Producător / HR Proizvođač / SL Proizvajalec / TR Üretici / MK Производител / SR Proizvođač / RU Производитель / VI Nhà sản xuất / ZH-CN 制造商 / ZH-TW 製造商
	EN Date of Manufacture / FR Date de fabrication / DE Herstellungsdatum / IT Data di produzione / ES Fecha de fabricación / PT Data de fabricação / NL Productiedatum / DA Fremstillingsdato / FI Valmistuspäivämäärä / SV Tillverkningsdatum / NO Produksjonsdato / EL Ημερομηνία κατασκευής / BG Дата на производство / ET Valmistamise kuupäev / LT Pagaminimo data / LV Ražošanas datums / PL Data produkcji / CS Datum výroby / SK Dátum výroby / HU Gyártás dátuma / RO Data fabricației / HR Datum proizvodnje / SL Datum proizvodnje / TR Üretim Tarihi / MK Датум на производство / SR Datum proizvodnje / RU Дата изготовления / ID Tanggal Produksi / VI Ngày sản xuất / ZH-CN 制造日期 / ZH-TW 製造日期
	EN Catalog number / FR Numéro de catalogue / DE Katalognummer / IT Numero di catalogo / ES Número de catálogo / PT Número de catálogo / NL Catalogusnummer / DA Katalognummer / FI Luettelonumero / SV Katalognummer / NO Katalognummer / EL Αριθμός καταλόγου / BG Каталоген номер / ET Katalooginumber / LT Katalogo numeris / LV Kataloga numurs / PL Numer katalogowy / CS Katalogové číslo / SK Katalógové číslo / HU Katalógusszám / RO Număr de catalog / HR Kataložki broj / SL Kataložka številka / TR Katalog numarasi / MK Број на каталог / SR Kataložki broj / RU Каталогный номер / ID Nomor katalog / VI Số catalog / ZH-CN 产品目录号 / ZH-TW 產品目錄編號
	EN Batch code / FR Numéro de lot / DE Chargennummer / IT Codice lotto / ES Código del lote / PT Código do lote / NL Partijcode / DA Batchkode / FI Eräkoodi / SV Batchnummer / NO Batchkode / EL Κωδικός παρτίδας / BG Партиден номер / ET Partii kood / LT Partijos kodas / LV Partijas kods / PL Kod partii / CS Kód šarže / SK Číslo šarže / HU Köteggkód / RO Codul lotului / HR Šifra serije / SL Serijska številka / TR Parti kodu / MK Сериски код / SR Šifra serije / RU Код партии / ID Kode batch / VI Mã lô / ZH-CN 批号 / ZH-TW 製造批號
	EN Caution, consult accompanying documents / FR Attention, consulter les documents joints / DE Achtung, Begleitdokumente lesen / IT Attenzione, consultare la documentazione allegata / ES Precaución, consultar los documentos adjuntos / PT Atenção, consultar os documentos em anexo / NL Let op, raadpleeg meegeleverde documenten / DA Forsigtighedsregler, se medfølgende dokumenter / FI Huomio, lue toimitetut asiakirjat / SV Viktigt! Läs medföljande dokument / NO Forsiktig, rådfør deg med vedlagte dokumenter / EL Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / BG Внимание, направте справка с придружаващата документация / ET Ettevaatust, lugege kaasnevaid dokumente / LT Dėmesio! Žr. pridedamus dokumentus / LV Uzmanību! Lasīt pavaddokumentus! / PL Uwaga: należy zapoznać się z załączonymi dokumentami / CS Pozor, seznamte se s přiloženými dokumenty / SK Upozornenie, pozrite si priložené dokumenty / HU Vigyázat, olvassa el a kísérfő dokumentumokat / RO Atenție, consultați documentele însoțitoare / HR Oprez, pročitati popratne dokumente / SL Previdno, preglejte priloženo dokumentacijo / TR Dikkat, birlikte verilen belgelere bakın / MK Внимание, консултирајте се со придружните документи / SR Oprez, pogledajte prateću dokumentaciju / RU Внимание! Прочтите сопроводительную документацию / ID Perhatian, lihat dokumen yang menyertai / VI Thận trọng, tham khảo tài liệu kèm theo / ZH-CN 注意, 请查阅随附文件 / ZH-TW 警告, 請查閱隨附文件說明
	EN Do not resterilize / FR Ne pas restériliser / DE Nicht erneut sterilisieren / IT Non risterrilizzare / ES No volver a esterilizar / PT Não reesterilizar / NL Niet opnieuw steriliseren / DA Må ikke resteriliseres / FI Ei saa steriloida uudestaan / SV Får inte resteriliseras / NO Må ikke resteriliseres / EL Μην επαναποστειρώνετε / BG Да не се стерилизира повторно / ET Mitte resteriliseerida / LT Pakartotinai nesterilizuoti / LV Nesterilizēt atkārtoti / PL Nie steryliżować ponownie / CS Nesterilizujte opakovaně / SK Nesterilizujte opakovane / HU Ne sterilizálja újra / RO A nu se resteriliza / HR Ne ponovno sterilizirati / SL Ponovna sterilizacija ni dovoljena / TR Tekrar sterilize etmeyin / MK He resterilizirajte / SR Ne sterilisati ponovo / RU He стерилизовать повторно / ID Jangan disterilisasi ulang / VI Không tái tiệt trùng / ZH-CN 请勿重复灭菌 / ZH-TW 請勿重新滅菌
	EN Do not re-use / FR Ne pas réutiliser / DE Nicht wiederverwenden / IT Non riutilizzare / ES No reutilizar / PT Não reutilizar / NL Niet hergebruiken / DA Må ikke genbruges / FI Ei saa käyttää uudelleen / SV Får inte återanvändas / NO Må ikke gjenbrukes / EL Μην επαναχρησιμοποιείτε / BG Да не се използва повторно / ET Mitte kordukasutada / LT Nenaudoti pakartotina / LV Nelietot atkārtoti / PL Nie używać ponownie / CS Nepoužívejte opakovaně / SK Na jednorazové použitie / HU Ne használja újra / RO A nu se reutiliza / HR Ne ponovno upotrebljavati / SL Ponovna uporaba ni dovoljena / TR Tekrar kullanmayın / MK He upotrebljavajte ponovno / SR Ne upotrebljavati više puta / RU He использовать повторно / ID Jangan gunakan ulang / VI Không tái sử dụng / ZH-CN 请勿重复使用 / ZH-TW 請勿重複使用

<p>STERILE EO</p>	<p>EN Sterilized using ethylene oxide / FR Stérilisé à l'oxyde d'éthylène / DE Sterilisiert mit Ethylenoxid / IT Sterilizzato con ossido di etilene / ES Esterilizado con óxido de etileno / PT Esterilizado com óxido de etileno / NL Gesteriliseerd met ethylenoxide / DA Steriliseret med ethylenoxid / FI Tämä tuote on steriloitu etyleenioksidilla / SV Steriliserad med etylenoxid / NO Sterilisert med etylenoksid / EL Αποστειρωμένο με τη χρήση αιθυλενοξειδίου / BG Стерилизирано с етиленов оксид / ET Steriliseeritud etüleenoksiidiga / LT Sterilizuota naudojant etileno oksidą / LX Sterilizēts ar etilēna oksīdu / PL Produkt sterylizowany tlenkiem etylenu / CS Sterilizováno etylenoxidem / SK Sterilizované pomocou etylénoxidu / HU Etilén-oxidall sterilizálva / RO Sterilizat cu oxid de etilenă / HR Sterilizirano etilen-oksidom / SL Sterilizirano z etilenoksidom / TR Etilen oksitle sterilize edilmiştir / MK Стерилизирано со употреба на етилен оксид / SR Sterilisan o etilen-oksidom / RU Стерилизовано этиленоксидом / ID Disterilkan dengan menggunakan etilena oksida / VI Sản phẩm đã được tiệt trùng bằng Ethylene Oxide / ZH-CN 经环氧乙烷灭菌 / ZH-TW 使用環氧乙烷滅菌</p>
	<p>EN Use-by date / FR Date limite d'utilisation / DE Verfallsdatum / IT Usare entro la data di scadenza / ES Fecha de caducidad / PT Prazo de validade / NL Uiterste gebruiksdatum / DA Anvendes inden / FI Viimeinen käyttöpäivä / SV Används senast / NO Utløpsdato / EL Ημερομηνία "Χρήση έως" / BG Срок на годност / ET Kõlblik kuni / LT Panaudoti iki nurodytos datos / LX Derīguma termiņš / PL Data ważności / CS Datum konce spotřeby / SK Spotrebujte do / HU Felhasználhatósági idő / RO Data expirării / HR Rok upotrebe / SL Datum izteka roka uporabnosti / TR Son kullanma tarihi / MK Датум на употреба / SR Rok upotrebe / RU Использовать до / ID Tanggal kedaluwarsa / VI Hạn sử dụng / ZH-CN 使用期限 / ZH-TW 使用期限</p> <p>EN Do not use this device after the indicated date (Year-Month-Day) / FR Ne pas utiliser ce dispositif après la date indiquée (année-mois-jour) / DE Nicht nach Ablauf des angegebenen Datums (Jahr-Monat-Tag) verwenden / IT Non usare questo dispositivo dopo la data indicata (anno-mese-giorno) / ES No utilizar este dispositivo después de la fecha indicada (Año-Mes-Día) / PT Não utilizar este dispositivo depois da data indicada (Ano-Mês-Dia) / NL Gebruik dit instrument niet na de aangegeven datum (jaar, maand, dag) / DA Anordningen må ikke anvendes efter den anførte dato (år-måned-dag) / FI Älä käytä laitetta ilmoitetun päivämäärän jälkeen (vuosi-kuukausi-päivä) / SV Använd inte denna enhet efter angivet datum (år-månad-dag) / NO Ikke bruk dette utstyret etter den angitte datoen (år-måned-dag) / EL Μη χρησιμοποιείτε αυτό το τεχνολογικό προϊόν μετά την ημερομηνία που επισμαίνεται (Έτος-Μήνας-Ημέρα) / BG Продуктът да не се използва след указаната дата (Година-месец-ден) / ET Ärge kasutage toodet pärast pakendil toodud kuupäeva (aasta-kuu-päev) / LT Šio prietaiso nenaudokite po nurodytos datos (metai-mėnuo-diena) / LX Neizmantojiet šo ierīci pēc norādītā datuma (Gads-Mēnesis-Datums) / PL Nie używać urządzenia po wskazanym dacie (Rok-Miesiąc-Dzień) / CS Tento prostředek nepoužívejte po uplynutí uvedeného data spotřeby (rok-měsíc-den) / SK Túto pomôcku nepoužívajte po uplynutí uvedeného dátumu (rok-mesiac-deň) / HU A feltüntetett dátum (év-hónap-nap) után az eszközt ne használja / RO Nu utilizați dispozitivul după data indicată (an-lună-zii) / HR Nemojte upotrebljavati ovaj proizvod nakon naznačena vijeka trajanja (godina-mjesec-dan) / SL Po izteku navedenega datuma pripomočka ne uporabljajte (leto-mesec-dan) / TR Bu cihaz belirlenen tarihten sonra kullanmayın (Yıl-Ay-Gün) / MK Не користете го овој уред по наведениот датум (година-месец-ден) / SR Ne koristiti ovo medicinsko sredstvo nakon označenog datuma (dan-mesec-godina) / RU Не использовать данное устройство позже указанной даты (год-месяц-день) / ID Jangan gunakan perangkat ini setelah melewati tanggal yang diindikasikan (Tahun-Bulan-Hari) / VI Không sử dụng thiết bị này sau ngày được chỉ định (Năm-Tháng-Ngày) / ZH-CN 超过指定日期 (年-月-日) 后请勿使用该器械 / ZH-TW 請勿在指定日期 (年-月-日) 後使用此裝置</p>
<p>MD</p>	<p>EN Medical Device / FR Dispositif médical / DE Medizinprodukt / IT Dispositivo medico / ES Dispositivo médico / PT Dispositivo médico / NL Medisch hulpmiddel / DA Medicinsk udstyr / FI Lääkinnällinen laite / SV Medicinteknisk produkt / NO Medisinsk utstyr / EL Ιατροτεχνολογικό προϊόν / BG Медицинско изделие / ET Meditsiiniseade / LT Medicinos priemonė / LX Mediciniskā ierīce / PL Wyrób medyczny / CS Zdravotnický prostředek / SK Zdravotnícka pomôcka / HU Orvostechnikai eszköz / RO Dispozitiv medical / HR Medicinski proizvod / SL Medicinski pripomoček / TR Tıbbi Cihaz / MK Медицински уред / SR Medicinsko sredstvo / RU Медицинское изделие / ID Alat Medis / VI Thiết bị Y tế / ZH-CN 医疗器械 / ZH-TW 醫療裝置</p>
<p>UDI</p>	<p>EN Unique Device Identifier / FR Identifiant unique du dispositif / DE Eindeutige Produktkennzeichnung (Unique Device Identifier) / IT Identificativo univoco del dispositivo / ES Identificador único del dispositivo / PT Identificador único de dispositivo / NL Unieke identificatiecode van het hulpmiddel / DA Unik enhedsidentifikationskode / FI Yksilöllinen laitetunniste / SV Unik produkt identifierare / NO Unik utstyrsidentifikator / EL Μοναδικός αναγνωριστικός κωδικός τεχνολογικού προϊόντος / BG Уникален идентификатор на изделието / ET Kordumatu identifitseerimiskood / LT Unikalusis priemonės identifikatorius / LX Unikālais ierīces identifikators / PL Unikalny identyfikator wyrobu / CS Jedinečný identifikátor prostředku / SK Unikátny identifikátor pomôcky / HU Egyedi eszközazonosító / RO Identificatorul unic al dispozitivului / HR Jedinstvena identifikacija proizvoda / SL Edinstveni identifikator pripomočka / TR Benzersiz Cihaz Tanımlayıcısı / MK Единствен идентификатор на уред / SR Jedinstveni identifikator medicinskog sredstva / RU Уникальный идентификатор изделия / ID Pengidentifikasi Unik Alat / VI Mã Định danh Thiết bị Duy nhất / ZH-CN 医疗器械唯一标识 / ZH-TW 唯一裝置識別碼</p>
	<p>EN Entity importing the medical device into the European Union / FR Entité qui importe le dispositif médical dans l'Union européenne / DE Unternehmen, das das Medizinprodukt in die Europäische Union einführt / IT Entità importatrice del dispositivo medico nell'Unione Europea / ES Entidad que importa el dispositivo médico a la Unión Europea / PT Entidade que importa o dispositivo médico para a União Europeia / NL Entiteit die het medische hulpmiddel in de Europese Unie importeert / DA Importør af det medicinske udstyr til den Europæiske Union / FI Yksikkö, joka tuo lääkitieteellisen laitteen Euroopan unioniin / SV Företaget som importerar den medicintekniska produkten till EU / NO Enhet som importerer det medisinske utstyret til EU / EL Οντότητα που εισάγει το ιατροτεχνολογικό προϊόν στην Ευρωπαϊκή Ένωση / BG Вносител на медицинското изделие в Европейския съюз / ET Üksus, kes impordib meditsiiniseadme Euroopa Liitu / LT Subjektas, importuojantis medicinos priemonę į Europos Sąjungą / LX Subjekts, kas importē medicīnisko ierīci uz Eiropas Savienību / PL Podmiot będący importerem wyrobu medycznego na terenie Unii Europejskiej / CS Subjekt dovážející zdravotnický prostředek do Evropské unie / SK Dovozca zdravotníckej pomôcky do Európskej únie / HU Az orvostechnikai eszköz az Európai Unióba importáló szervezet / RO Entitatea care importă dispozitivul medical în Uniunea Europeană / HR Subjekt koji uvozi medicinski proizvod u Europsku uniju / SL Subjekt, ki medicinski pripomoček uvaža v Evropsko unijo / TR Tıbbi cihaz Avrupa Birliği'ne ithal eden kuruluş / MK Субјект кој го увезува медицинскиот уред во Европската Унија / SR Označava pravnog subjekta koji uvozi medicinsko sredstvo u Evropsku uniju / RU Организация-импортер медицинского изделия в Европейский союз / ID Entitas yang mengimport alat medis ke Uni Eropa / VI Tổ chức nhập khẩu thiết bị y tế vào Liên minh châu Âu / ZH-CN 进口医疗器械到欧盟的实体 / ZH-TW 將醫療裝置進口到的實體歐盟</p>

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Coronary Microcatheter

Table of Contents

1. DEVICE DESCRIPTION
 - 1.1 Device Components
2. INTENDED PURPOSE
3. INDICATIONS
4. CONTRAINDICATIONS
5. INTENDED USERS
6. USE IN SPECIAL POPULATIONS
7. INTENDED CLINICAL BENEFIT
8. WARNINGS
9. PRECAUTIONS
10. OPERATOR'S MANUAL
 - 10.1 Materials Required
 - 10.2 Inspection Prior to Use
 - 10.3 Instructions for Use
11. POTENTIAL ADVERSE EVENTS
12. PATIENT COUNSELING INFORMATION
13. HOW SUPPLIED AND DISPOSAL
14. SUMMARY REPORT ON SAFETY AND CLINICAL PERFORMANCE
15. WARRANTY

1. DEVICE DESCRIPTION

Biosensors BioMC™ Coronary Microcatheter is an Over-the-Wire microcatheter that is designed to provide guidewire support during lesion crossing, facilitate guidewire exchange and deliver contrast media. It is available in usable catheter lengths of 135 and 150 cm. The microcatheter is packed in a holder tube and sealed within a sterile pack.

1.1. Device Components

- The microcatheter consists of a distal tip, a shaft tube, a strain relief and a hub.
- The distal tip and the distal part of the shaft are coated with a lubricious hydrophilic solution.
 - The inner lumen is made using a fluoropolymer for smooth guidewire exchange.
 - The distal end of the microcatheter has a radiopaque tip to distinguish the end.

Table 1: Device Description

Microcatheter Working Length	135 cm	150 cm
Inner Liner Material	Polytetrafluoroethylene (PTFE)	
Microcatheter Shaft Profile	Proximal: 2.8 F / 0.037" / 0.94 mm Distal: 2.3 F / 0.031" / 0.78 mm	
Maximum Injection Pressure	2070 kPa (300 PSI)	
Guiding Catheter Compatibility	5F	
Guidewire Inner Lumen Compatibility	0.014"	

Table 2: Product range available

Product Code	Working Length (cm)
MC1-135	135
MC1-150	150

2. INTENDED PURPOSE

The BioMC Coronary Microcatheter is intended to be used on patients eligible for Percutaneous Coronary Intervention (PCI) to treat Coronary Arterial Disease (CAD).

3. INDICATIONS

The BioMC Coronary Microcatheter is indicated for use:

- in the coronary vasculature to provide support to facilitate both the placement and exchange of guidewires and to deliver contrast media in patients eligible for Percutaneous Coronary Intervention (PCI).
- in patients eligible for PCI of Chronic Total Occlusions (CTOs) to provide additional guide wire support and increase guide wire penetration force to facilitate crossing of the occluded segment.

4. CONTRAINDICATIONS

- Patients with known sensitivity to contrast agents that cannot be treated prophylactically prior to the use of the microcatheter.

5. INTENDED USERS

The microcatheter is a professional use medical device that is intended to be used by professionals who have received appropriate training and education in Percutaneous Coronary Intervention (PCI) procedures. The microcatheter should only be used by interventional cardiologists. There is no specific procedural or training requirements for this product to be provided by Biosensors to the healthcare professionals.

6. USE IN SPECIAL POPULATIONS

The safety and effectiveness of the microcatheter has not been established in the following patient populations:

- Pregnancy: there is no data available for use in pregnant women.
 - Pediatric use: the safety and efficacy has not been established.
- Carefully consider whether it is appropriate to use the microcatheter in the above patient populations.

7. INTENDED CLINICAL BENEFIT

The clinical benefits of the microcatheter are expected to be the same as for the other coronary microcatheters.

- The microcatheter is expected to enhance the deliverability and penetration force of coronary guide wires in patients undergoing PCI, and thus facilitate lesion crossing.
- Myocardial revascularization by PCI is a minimally invasive procedure that improves blood flow to the heart, decreasing heart-related chest pain (angina) and improving

exercise capacity¹.

- In Acute Coronary Syndromes, PCI reduces the risk of death and recurrent Myocardial Infarction².
- In Chronic Total Occlusion (CTO), PCI improves angina symptoms and Quality of Life (QoL)³.

8. WARNINGS

- Ensure that the packaging has not been damaged or opened as this may indicate a breach of the sterile barrier. Check the "Use by date" indicated next to the representative symbol on the labels and do not use a product that has reached or exceeded its labeled expiration date.
- The sterility and stability of the device cannot be guaranteed once the pouch has been opened and hence the device MUST be used promptly.
- Un-used devices should be returned to Biosensors International™⁴ and should not be re-stocked.
- This device is intended for single use only and must not be reused in another procedure. Do Not resterilize the device as this can compromise its performance, can lead to device failure and procedure complications with severe injury or patient death. Reuse, reprocessing and reesterilization bear the risk of cross contamination and patient to patient infection.
- Inappropriate storage of the device may compromise its performance, can lead to device failure and procedure complications with severe injury or patient death.
- Only physicians who have received appropriate training and education should perform PCI procedures and the risk/benefit of each procedure shall be assessed for each patient by the physician.
- PCI should only be performed in healthcare settings that allow emergency coronary artery bypass graft surgery (CABG) to be performed in the event of a potentially life-threatening or other serious complication.
- The microcatheter must not be used for drug infusion, except for the contrast media. The microcatheter is not designed for drug infusion and its safety has not been established.
- Do not modify this microcatheter for any reason. Use of a modified product may cause damage to blood vessels and/or accidents.
- Excess rotational force should not be applied if the microcatheter is bent. The microcatheter may be damaged.
- Do not insert the guide wire by force or advance it rapidly when the microcatheter is bent or twisted. Such manipulations may cause breakage or damage of the microcatheter, or perforation of the blood vessel.
- Advancement of the microcatheter must always be performed using high-resolution fluoroscopic guidance. Particular attention should be paid when inserting or withdrawing the microcatheter across stenotic and/or occluded segments particularly in vessels with a diameter less than the product. The use of excess force can result in damage or separation of the microcatheter and may cause vascular injury and/or perforation, possibly leading to a life-threatening adverse event.
- If any abnormal resistance is felt during use of the microcatheter, do not continue the procedure until the cause is established. If it is suspected that the microcatheter is not functioning correctly, avoid excessive manipulations, and carefully remove the entire microcatheter system while paying full attention to avoid complications. Continuing the procedure while the cause of the problem has not been identified may cause damage to or separation of the microcatheter, damage to the blood vessel and life-threatening adverse events.
- Always advance the guide wire ahead of the microcatheter before attempting any manipulation of the microcatheter. If the guide wire is not advanced ahead of the microcatheter, the blood vessel may be damaged or perforated, or the microcatheter may be damaged.
- Always hold the connector with one hand and turn the microcatheter carefully while regularly releasing the accumulated torsion of the microcatheter. Never turn the microcatheter continuously while holding the connector with both hands and never use any other means to apply force. When releasing the accumulated torsion, be sure to open the hemostatic valve on the Y-connector. Do not turn the microcatheter in the same direction, either clockwise or counterclockwise, for more than 15 consecutive turns. Continuing rotation may damage or break the microcatheter or damage the blood vessels leading to life threatening adverse events. If resistance is felt while turning the microcatheter, do not proceed with further rotation even if the 15-turn limit has not been reached. Identify the cause

1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.

2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.

3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.

4 Please contact the Sales and Customer Service of your region or local distributor for return of goods.

of resistance under fluoroscopy and take appropriate action. Never continue the operation without identifying the cause.

- This microcatheter is coated with hydrophilic coating and is therefore highly lubricious. Always confirm the position of the distal end of this microcatheter and manipulate it carefully to avoid uncontrolled forward movement.
- Do not infuse contrast media when the microcatheter is bent or occluded as this may cause over-expansion and breakage of the device.
- When injecting contrast media, the injection pressure must not exceed 2070 kPa (300psi). Exceeding the maximum injection pressure may cause damage to the microcatheter.
- When injecting contrast media, use fluoroscopy to confirm that the contrast media is being infused from the tip of the device. If no contrast media is visualized the infusion must be stopped and the device must be replaced with a new one. If the device lumen is occluded, it may be damaged, over-expanded or ruptured even at less than the maximum injection pressure, resulting in a life-threatening adverse event.
- If the microcatheter is being used for guide wire exchange, exercise care when removing the original guide wire and inserting the new one. If any abnormal resistance is felt, the procedure should be discontinued immediately, and the microcatheter and the guide wire should be withdrawn together.
- Repeated insertion and withdrawal of the device may lead to deterioration of the hydrophilic coating. Continuous use of a device with a deteriorated hydrophilic coating may cause vascular damage and may increase the risk of the microcatheter becoming trapped, resulting in a life-threatening adverse event due to damage and/or separation of tip.

9. PRECAUTIONS

- When inserting the guide wire into a microcatheter which is already placed in the blood vessel, carefully advance the guide wire to minimize the risk of damage to the microcatheter within tortuous or angulated segments.
- Confirm that the microcatheter does not have a kink, knot, torsion, or occlusion before injecting contrast media.
- Select the appropriate size of guiding catheter and guide wire to use in combination with this microcatheter. Refer to product labelling.
- When using a guiding catheter fitted with a stopcock, do not manipulate the stopcock after inserting the microcatheter into the guiding catheter. If the stopcock is manipulated during insertion, the microcatheter may be damaged.
- Manipulate the microcatheter carefully to avoid damage, kinking, or bending, especially when inserting into the guiding catheter.
- Check the patient's condition before the procedure. Provide appropriate anticoagulant therapy as per routine standards of care for patients undergoing PCI.
- The surface of this microcatheter is coated with hydrophilic polymer. Flush the surface and the lumen of the microcatheter continuously with heparinized and sterilized saline during its use to maintain lubricity. Flush the lumen of the microcatheter sufficiently with heparinized and sterilized saline especially after injecting contrast media.
- When using a hemostatic adaptor (i.e., Y-connector), avoid excessive tightening of the hemostatic valve as the microcatheter may be damaged.
- Since this microcatheter is not designed for use in combination with a power injector, do not use it in combination with a power injector.
- Do not expose the device to organic solvents, e.g. isopropyl alcohol. Such an exposure can degrade the microcatheter performance.

10. OPERATOR'S MANUAL

Any serious incident in relation to the device should be reported to Biosensors International™ and the relevant competent authority as required by local regulations.

10.1. Materials Required

One or more of each of the following material are required but not supplied with the microcatheter:

- Arterial sheath and dilatator set
- Guiding catheter (femoral or brachial) in the appropriate configuration to selectively cannulate the coronary artery and with minimum inner diameter of 0.056"/1.42 mm (5F compatible)
- Guidewire, 0.014" / 0.36 mm maximum diameter x 190 cm minimum length
- Guidewire torque device
- Hemostatic adaptor
- Three-way stopcock
- Luer-lock syringes
- Sterile saline or heparinized sterile saline
- Contrast medium diluted 1:1 with normal saline

10.2. Inspection Prior to Use

Prior to use, inspect the sterile package integrity, do not use if packaging is damaged. Carefully examine the microcatheter prior to use to detect any bends, kinks, or other damage. Verify all equipment as well to be used during the procedure. Do not use any damaged equipment, a product from a damaged package or breached sterile pouch or an expired product. Verify that the microcatheter length is suitable for the specific procedure for which it is intended.

10.3. Instructions for Use

Prepare each equipment to be used following instructions of its manufacturer.

Complete the following steps to use the microcatheter:

1. Using aseptic technique, carefully open the pouch and remove the holder tube containing the microcatheter from the sterile packaging.
2. Remove the microcatheter from the holder tube and inspect it to verify that the microcatheter is free of damage prior to insertion.
3. Attach a syringe filled with heparinized saline to the connector of the microcatheter and inject enough fluid to purge the air from the inside of the microcatheter.
4. Wet the microcatheter surface entirely with heparinized saline solution or sterile water to activate the hydrophilic coating.
5. Insert a compatible guide wire through the connector and bring the tip of the guide wire in line with the tip of this microcatheter. If the guide wire is inserted through the tip of this microcatheter, care should be taken not to cause any damage to the microcatheter. Also, if the microcatheter is bent or kinked, discontinue its use. If the microcatheter is kinked it may cause severe damage to the patient.
6. Introduce the microcatheter and guide wire assembly into the guiding catheter via the hemostatic adaptor (if used) and advance the assembly until it reaches 2 to 3 cm proximal of the tip of the guiding catheter.
7. Using fluoroscopy, advance this microcatheter and guide wire assembly into the vascular system until it is close to the stenotic area, making sure the guide wire is always ahead of the microcatheter. Advance the guide wire carefully until it passes the target area. Continue advancing the guide wire as distal as possible into the blood vessel, and once it is placed there, check the position by imaging from the guiding catheter. The position of the guide wire must be checked by imaging from multiple angles to confirm that the guide wire is definitely inserted into the target blood vessel.
8. After loosening the hemostatic valve, hold the guide wire and guiding catheter firmly. Then advance this microcatheter gradually along the guide wire until the tip has passed through the stenotic area, using the tip of this microcatheter as a guide. Procedures inside the blood vessel should be conducted with care, because this microcatheter is hydrophilic coated.
9. The user can rotate the device when inserting, withdrawing, and passing through stenotic areas.
CAUTION: Do not turn the microcatheter in the same direction, either clockwise or counterclockwise, for more than 15 consecutive turns. If the device is trapped or suspected to be trapped, rotating operation must be avoided.
10. For guide wire exchange, withdraw the initial guide wire from microcatheter. Then, insert the new guide wire via the hemostatic adaptor (i.e., Y-connector). Using fluoroscopic guidance, advance the new guide wire to the end-hole of microcatheter.
11. In case injection of contrast media is required, withdraw the guide wire from the microcatheter and inject the contrast media from the microcatheter hub with a small volume luer lock syringe.
CAUTION: slowly inject a small amount of contrast media and observe the outflow from the distal end of the product to ensure the lumen is patent. If no contrast is seen exiting the microcatheter, do not continue the injection as this may result in deformation or damage to the product.
12. When removing the microcatheter, loosen the hemostatic valve of the hemostatic adaptor. Remove this microcatheter while keeping the guide wire stable in the blood vessel. When this microcatheter is removed, check the position of the guide wire under fluoroscopy. Also, if any resistance is felt during the removal of this microcatheter, remove all devices including the parent microcatheter and the guide wire. After removal of this microcatheter, tighten the hemostatic valve of the hemostatic adaptor.
13. After completing the procedure, withdraw the microcatheter immediately and discard.

11. POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a microcatheter in native coronary arteries, include but not limited to:

- Cardiac complications: acute myocardial infarction, ischemia, heart failure
- Coronary artery complications: abrupt vessel closure or spasm, aneurysm, pseudoaneurysm, dissection, perforation, rupture, injury of the treated coronary

artery, vascular occlusion

- Arrhythmic complications: arrhythmias
- Respiratory complications: embolism
- Vascular access complications: arteriovenous fistula, hemorrhage, infection
- Neurologic complications: stroke, Transient Ischemic Attack (TIA)
- Allergic reaction to anti-coagulation and/or anti-thrombotic therapy, contrast medium, or delivery system materials
- Death
- Hypotension, hypertension
- Renal failure

12. PATIENT COUNSELING INFORMATION

Physicians should provide the patient with information on risks/benefits (e.g., contraindications, warnings, precautions, adverse events and complications and alternative treatments), follow-up visits and medication to be taken.

13. HOW SUPPLIED AND DISPOSAL

STERILE, NON-PYROGENIC: Package contents are sterile unless package is opened or damaged. Do not use if package is open or damaged. This device is sterilized with ethylene oxide sterilization and is non-pyrogenic.

CONTENTS: One (1) Microcatheter and one (1) Instruction for Use.

STORAGE: Store in a cool, dark and dry place. Do not store above 25°C.

DISPOSAL: Disposal of the device and its packaging should be handled in an environmentally sustainable manner according to local regulations. Hazardous waste from contaminated devices and packaging may present a biohazard and must be disposed of in appropriate containers which meet specific technical requirements.

14. SUMMARY REPORT ON SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the BioMC Coronary Microcatheter can be found in the European database on medical devices (Eudamed): <https://ec.europa.eu/tools/eudamed>, by entering the Basic UDI-DI of the microcatheter (8888893MC1RR).

15. WARRANTY

Biosensors International and their respective affiliates (collectively Biosensors) warrants that its products are manufactured to the specifications set forth on its packaging, instructions for use and related literature.

This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied, by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.

Biosensors neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product.

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² Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020; 13(2):e006363.

³ Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr; 36(2):145-157.

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- 1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 - 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020; 13(2):e006363.
 - 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr; 36(2):145-157.
- 4 Entre em contacto com o serviço de vendas e apoio ao cliente da sua região ou com o seu distribuidor local para a devolução de produtos.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.
 4 Neem contact op met de Verkoop- en Klantenservice in uw regio of lokale distributeur voor de retourzending van goederen.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.
 4 Kontakta försäljnings- och kundservice för din region eller din lokala distributor beträffande retur av varor.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.
 4 Επικοινωνήστε με το Τμήμα Πωλήσεων και Εξυπηρέτησης Πελατών της περιοχής σας ή με τον τοπικό διανομέα για την επιστροφή των προϊόντων.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.
 4 За връщане на стоките, се свържете с отдел „Продажи и обслужване на клиенти“ във вашия регион или с местния дистрибутор.

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2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. *Circ Cardiovasc Qual Outcomes*. 2020; 13(2):e006363.

3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. *Cardiovasc Interv Ther*. 2021 Apr;36(2):145-157.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020; 13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr; 36(2):145-157.
 4 W celu zwrotu produktu należy skontaktować się z regionalnym działem sprzedaży i obsługi klienta lub lokalnym dystrybutorem.

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2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020; 13(2):e006363.

3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr; 36(2):145-157.

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² Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.

³ Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.

⁴ Za povrat proizvoda obratite se službi za prodaju i korisničku podršku u svojoj regiji ili lokalnom distributeru.

mikrokatereta, odmah prestanite okretati mikrokateret čak i ako nije dosegnuto ograničenje od 15 uzastopnih okretaja. Utvrdite uzrok otpora s pomoću fluoroskopije i poduzmite odgovarajuće mjere. Nemojte nastaviti okretati mikrokateret ako niste utvrdili uzrok otpora.

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- 1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
- 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020;13(2):e006363.
- 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr;36(2):145-157.
- 4 контактирајте со службата за продажба и клиенти на вашиот регион или локалниот дистрибутер за враќање на стоката.

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1 Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J. 2019 Jan 7; 40(2):87-165.
 2 Chacko et al. Effects of Percutaneous Coronary Intervention on Death and Myocardial Infarction Stratified by Stable and Unstable Coronary Artery Disease: A Meta-Analysis of Randomized Controlled Trials. Circ Cardiovasc Qual Outcomes. 2020; 13(2):e006363.
 3 Morino et al. A contemporary review of clinical significances of percutaneous coronary intervention for chronic total occlusions, with some Japanese insights. Cardiovasc Interv Ther. 2021 Apr; 36(2):145-157.
 4 Для возврата изделий обратитесь в службу продаж и обслуживания клиентов в вашем регионе или свяжитесь с местным дистрибутором.

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THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

BIOSENSORS
INTERNATIONAL™

Biosensors Europe SA
Rue de Lausanne 29
1110 Morges
Switzerland
Tel: +41 21 804 8000
Fax: +41 21 804 8001
www.biosensors.com

12248-000 rev02
2026-01

CE
0344

BioMC is a trademark or registered trademark
of Biosensors International Group, Ltd.

© 2026 - Biosensors International Group, Ltd. All rights reserved

BioMC
EUROPEAN MEGACALISTER