

BIPOLAR

Pacing Catheter



**BIOSENSORS
INTERNATIONAL**

Improving Lives Through Medical Innovation

Uncompromised Quality and Reliability

Drawing upon our sound track record as a specialised pulmonary artery catheter manufacturer, Biosensors International is complementing the existing range of right heart catheters with a new addition - bipolar pacing catheter.

The development of the bipolar pacing catheter for temporary transvenous cardiac pacing brings to bear Biosensors's experience in this field.


The launch of this product is in response to market demand for a safe, secure and accurate bipolar pacing catheter to support critical procedures in the Cath Lab, CCU and Emergency Room settings.

Safe, Sure and Accurate

When safety and accuracy are paramount during critical procedures, cardiologists and electro-physiologists know they can count on Biosensor's bipolar pacing catheters.

Medical practitioners can work with confidence in determining the depth of catheter insertion as the length of the tubing is marked in 10cm increments. Radiopaque polyurethane tubing facilitates the use of fluoroscope for greater accuracy in catheter insertion and in verifying the position of the catheter.

Added safety is ensured by a lumen incorporated into the catheter to provide electrical conductor isolation from electrodes to the connector pins. An optional inflation balloon - located at the distal tip - can be inflated to enable safe flotation of the catheter to the right ventricle.



Specifications

Product Code	Catheter Size	Catheter Length	Inflated Balloon Capacity	Maximum Inflation Diamete	Electrode Width	Electrode Spacing	Contamination Shield	Recommended Introducer Size
With Balloon								
BP2502-10	5F	110cm	1cc	9mm	3mm	10mm	NA	6F
BP2502CS-10	5F	110cm	1cc	9mm	3mm	10mm	110cm	6F
Without Balloon								
BPX2502-10	5F	110cm	NA	NA	3mm	10mm	NA	6F
BPX2502CS-10	5F	110cm	NA	NA	3mm	10mm	110cm	6F

CE0044

Singapore:
Biosensors International Pte Ltd
 21 Kallang Avenue, #07-165/171,
 Singapore 339412
 Tel: +65 6411 0500
 Fax: +65 6298 6242
 Email: sales@biosensors.com.sg

United States of America:
Biosensors International USA
 20250 Acacia Street,
 Suite 200, Newport Beach,
 CA 92660, USA
 Tel: +1 (949) 553 8300
 Fax: +1 (949) 553 9129
 Email: sales@biosensors-usa.com



www.biosensors.com

Netherlands:
Biosensors BV
 Weeresteinstraat 147,
 2182 GV, Hillegom, The Netherlands
 Tel: +31 (25) 251 7676
 Fax: +31 (25) 252 6782
 Email: sales-ccp-bv@biosensors.com

Japan:
Biosensors Japan Co., Ltd.
 Bellebs Nagayama 501,
 1-5 Nagayama, Tama-Shi,
 Tokyo 206 0025, Japan
 Tel: +81 (42) 355 3311
 Fax: +81 (42) 355 3312
 Email: cs@sj.co.jp

China:
Biosensors International China Sales Office
 Room 408, Guobin Building,
 No. 11, Fuchengmenwai Street, Xicheng District,
 Beijing 100037, China
 Tel: +86 (10) 6800 3481
 Fax: +86 (10) 6800 4202
 Email: sales.icu@jwmsgrp.com
 chifei.zhou@jwmsgrp.com