



## Embolectomy Catheter

- Biosensors Embolectomy Catheters offer simple, rapid removal of emboli from the arterial system.
- Proprietary "low-profile" winding technique for balloon attachment insures balloon symmetry, providing uniform contact with vessel walls.
- Soft, rounded tip promotes easy arterial insertion for reduced trauma, arteriopuncture and plaque avulsion.
- Flexible polyvinylchloride body with "workhorse" latex balloon.
- Winged hub for easy catheter manipulation during insertion and attachment of balloon inflation syringe.
- Available in 4 lengths and 6 diameters, ranging from 2 French to 7 French.
- Size and balloon inflation capacity printed on proximal end for convenient and rapid identification.

### Specifications

Product Code	Catheter Size	Length (cm)	Colour Code	Maximum Liquid Volume (cc)	Maximum CO <sub>2</sub> Volume (cc)	Inflated Balloon Diameter (mm)	Maximum Pull Force Inflated Balloon (lb)
Specifications							
EL 3260	2F	60	Purple	0.08	0.2	4	0.5
EL 3280	2F	80	Purple	0.08	0.2	4	0.5
EL 3340	3F	40	Green	0.20	0.6	6	0.7
EL 3380	3F	80	Green	0.20	0.6	6	0.7
EL 1440	4F	40	Red	0.75	1.7	9	1.5
EL 1480	4F	80	Red	0.75	1.7	9	1.5
EL 14100	4F	100	Red	0.75	1.7	9	1.5
EL 1540	5F	40	White	1.50	3.0	11	2.0
EL 1580	5F	80	White	1.50	3.0	11	2.0
EL 1640	6F	40	Blue	2.00	4.5	13	2.5
EL 1680	6F	80	Blue	2.00	4.5	13	2.5
EL 1740	7F	40	Yellow	2.50	5.0	14	3.5
EL 1780	7F	80	Yellow	2.50	5.0	14	3.5

Note: Only 2F and 3F Embolectomy Catheters are nylon body

CE 0044

Singapore:  
**Biosensors International Pte Ltd**  
 21 Kallang Avenue, #07-165/171,  
 Singapore 339412  
 Tel: +65 6411 0500  
 Fax: +65 6298 6242  
 Email: sales@biosensors.com.sg



[www.biosensors.com](http://www.biosensors.com)

United States of America:  
**Biosensors International USA**  
 20250 Acacia Street,  
 Suite 200, Newport Beach,  
 CA 92660, USA  
 Tel: +1 (949) 553 8300  
 Fax: +1 (949) 553 9129  
 Email: sales@biosensors-usa.com

China:  
**Biosensors International China Sales Office**  
 Room 408, Guobin Building,  
 No. 11, Fuchengmenwai Street, Xicheng District,  
 Beijing 100037, China  
 Tel: +86 (10) 6800 3481  
 Fax: +86 (10) 6800 4202  
 Email: sales.icu@jwmsgrp.com

Japan:  
**Biosensors Japan Co., Ltd.**  
 Bellebs Nagayama 501,  
 1-5 Nagayama, Tama-Shi,  
 Tokyo 206 0025, Japan  
 Tel: +81 (42) 355 3311  
 Fax: +81 (42) 355 3312  
 Email: cs@sj.co.jp

Netherlands:  
**Biosensors BV**  
 Weeresteinstraat 147,  
 2182 GV, Hillegom, The Netherlands  
 Tel: +31 (25) 251 7676  
 Fax: +31 (25) 252 6782  
 Email: sales-ccp-bv@biosensors.com

## Embolectomy Catheter



**BIOSENSORS  
INTERNATIONAL**

Improving Lives Through Medical Innovation